

Smart Transmitter/Gas Detector Head

SD-10X (TYPE H)

Operating Manual

(PT2E-177)

(PT2-177)

RIKEN KEIKI Co., Ltd.

2-7-6 Azusawa, Itabashi-ku, Tokyo, 174-8744, Japan

Phone: +81-3-3966-1113 Fax: +81-3-3558-9110

E-mail: intdept@rikenkeiki.co.jp Web site: https://www.rikenkeiki.co.jp/

Operating Precautions

This detector is an oxygen deficiency detector that detects oxygen in the air and triggers an oxygen deficiency alarm.

The oxygen deficiency detector is a safety unit, not an analyzer or densitometer which performs quantitative/qualitative analysis/measurement for oxygen.

Please fully understand the following points before using it, so that it can be used properly.

- 1. The readings of the detector fluctuate slightly in response to changes in the air pressure. In particular, be careful of alarm activation when a low air pressure is brought in by typhoon. In addition, it may be fluctuated by environmental (temperature, humidity etc.) changes in the installation site.
- 2. The alarm must be set within a range where the performance of the detector can be ensured. The standard alarm setpoint value is 18 vol%.
- 3. If use in safety unit, the analog signal output of the detector must be used for an indicator or external recorder.

If these outputs are used to control other units, we shall not be responsible for any malfunctions.

4. For maintenance of the detector, it must go through a regular maintenance, including replacement and adjustment of the regular replacement parts as specified in the operating manual. In addition, because this is a safety unit, it is recommended that a regular maintenance and span adjustment are performed every six months.

<Contents>

1 Outline of the Product	2 2 2
2 Important Notices on Safety 2-1. Danger cases 2-2. Warning cases 2-3. Precautions 2-4. Safety Information	4 5 6
3 Product Components	9 10
4 How to Use	12 13 14 16
5 How to Operate 5-1. Preparation for start-up 5-2. Basic operating procedures 5-3. How to start the gas detector 5-4. Modes 5-5. User mode 5-6. How to exit	21 21 22 23
6 Operations and Functions 6-1. Alarm activation 6-2. Fault alarm activation 6-3. External output operation 6-4. Other functions	27 28 29
7 Maintenance	31 33 40
8 Storage, Relocation and Disposal	45 45
9 Troubleshooting	46
10 Product Specifications 10-1. List of specifications 10-2. List of accessories 10-3. Detection principle	48 51
11 Definition of Terms	53

1 Outline of the Product 1-1. Preface

1

Outline of the Product

1-1. Preface

Thank you for choosing our smart transmitter/gas detector head SD-1OX(TYPE H). Please check that the model number of the product you purchased is included in the specifications on this manual.

This manual explains how to use the detector and its specifications. It contains information required for using the gas detector properly. Not only the first-time users but also the users who have already used the product must read and understand the operating manual to enhance the knowledge and experience before using the detector.

Product specifications may be abbreviated in this document as follows.

Japan Ex specification : Japan specification

ATEX / IECEx / UKEX specification : Export specification

1-2. Purpose of use


- This detector is a fixed type gas detector head that detects oxygen and performs the alarm activation when the gas concentration is over the setting value.
- The detector detects abnormalities in the air caused by presence of gases or other reasons (leak) with the built-in gas sensor. The concentrations of detected gases are displayed on the seven-segment LED.
- The detector outputs gas concentration in 4 20 mA.
- The detector has HART communication function.

1-3. Definition of DANGER, WARNING, CAUTION and NOTE


DANGER	This message indicates that improper handling may cause serious damage on life, health or assets.
WARNING	This message indicates that improper handling may cause serious damage on health or assets.
This message indicates that improper handling may cause minor of health or assets.	
NOTE	This message indicates advice on handling.

1-4. Method of confirmation for Standards and Explosion proof specification

This instrument has some specification depends on standard and explosion proof certificate. Please confirm the detector specification before using. Please refer Declaration of Conformity that is at the end of this manual if you have CE/UKCA marking type. You can confirm instrument specification to see name plate as follows.



ATEX/IECEx/UKEX, CE/UKCA marking type name plate



Japan Ex type name plate

2

Important Notices on Safety

2-1. Danger cases



DANGER

About explosion-proof

- The window plate material is a polycarbonate resin. Do not use organic solvents and alkali
 types (liquid or vapor). It may cause the color and shape of the window plate to be changed.
- Do not open the lid when applying current.
- Do not attempt to repair the detector.
- For the lid, use hexagon socket head bolts specified by RIKEN KEIKI.
- Do not apply a strong force or shock to the window plate. The explosion-proof performances may be deteriorated due to damages. The explosion protection condition for the detector is "low (2J)" possibility of mechanical damage.
- Ambient temperature:

Export Specification: -20 to +60°C Japan Specification: -20 to +53°C

(The ambient temperature on explosion-proof certification is described above and it indicates the temperature range which can maintain the explosion-proof performance and not the product performance, see "10-1. List of specifications" for details.)

- Do not replace parts at your sole discretion but contact RIKEN KEIKI if the transparent window
 has a crack or the explosion-proof joint surface is abnormal, or the clamping screw or bolt is
 changed, lost etc.
- Do not repair the explosion-proof joint surface.
- Be careful not to hit sharp tools against the window plate.

2-2. Warning cases



WARNING

Power supply

Before turning on the detector, always check that the voltage is properly applied. Do not use an unstable power supply because it may cause malfunctions.

Need of grounding circuit

Do not cut the grounding circuit or disconnect the wire from the grounding terminal.

Defects in protective functions

Before starting the detector, check the protective functions for defects. When seeming defects are found in the protective functions, such as protective grounding, do not start the detector.

External connection

Before connecting the detector to the external control circuit, securely connect it to a protective grounding circuit.

Span adjustment in atmosphere

When the span adjustment is performed in the atmosphere, check the atmosphere for freshness before beginning the adjustment. If interference gases exist, the adjustment cannot be performed properly, thus leading to dangers when the oxygen concentration fluctuates.

Response to gas alarm

Issuance of a gas alarm indicates that there are extreme dangers. Take proper actions based on your judgment.

Handling of sensor

Do not disassemble the sensor used by the detector because it contains electrolyte. Electrolyte may cause severe skin burns if it contacts skin, while it may cause blindness if its contacts eyes. If electrolyte is adhered on your clothes, that part on your clothes is discolored or its material is decomposed. If contact occurs, rinse the area immediately with a large quantity of water.

2-3. Precautions



CAUTION

Do not use a transceiver near the detector.

Radio wave from a transceiver or other radio wave transmitting device near the detector or its cables may disturb readings. If a transceiver or other radio wave transmitting device is used, it must be used in a place where it disturbs nothing.

To restart the detector, wait for five seconds or more before doing it.

Restarting the detector in less than five seconds may cause errors.

Do not use the external output of the detector to control other units.

This is not a control unit. It is not allowed to use the external output of the detector to control other units.

Do not disassemble/modify the detector, or change the settings if not necessary.

Disassembling/modifying the detector will invalidate the warranty of the performance. Changing the settings without understanding the specifications may cause alarm malfunctions. Please use the detector properly in accordance with the operating manual.

Avoid applying organic solvents and others to the window plate for a long time.

The window plate material is a polycarbonate resin. When organic solvents (liquid or highly-concentrated vapor) and others are applied to the plate for a long time, its color and shape may be changed.

Never fail to perform a regular maintenance.

Since this is a safety unit, a regular maintenance must be performed to ensure safety. Continuing to use the detector without performing a maintenance will deteriorate the sensitivity of the sensor, thus resulting in inaccurate gas detection.

2-4. Safety Information

Necessary information for explosion proof construction of Model SD-1OX.

The Model SD-1OX is a fixed mount, continuous-monitoring detector head and provides a 4-20mA signal which indicates the target gas reading for use by a gas monitoring controller, recording device, or programmable controller.

DEKRA21UKEX0368X

<ATEX/IECEx/UKEX Specification>

Technical Data

(Protection Method)Flameproof enclosure "d"(Certificate Numbers)IECEx DEK 13.0024XDEKRA 13ATEX0035 X

(Group) II (Category) 2G

(Type of Protection and Marking code) Ex db IIC T6

(Equipment Protection Level) Gb

(Ambient Temperature)* -20°C to +60°C

(Electrical Data) Supply voltage: DC24 V (Typ.)


Output signal (4 to 20 mA): DC24V 22mA

(Applicable Standard) IEC 60079-0: 2017, IEC 60079-1: 2014-06 EN IEC 60079-0: 2018, EN 60079-1: 2014

BS EN IEC 60079-0:2018, BS EN 60079-1:2014

(Manual Number) PT2E-177

Installation



^{*} The ambient temperature on explosion-proof certification is described above and it indicates the temperature range which can maintain the explosion-proof performance and not the product performance, see "10. Product specifications" for details.

<Japan Ex Specification>

Technical Data

(Protection Method) Flameproof enclosure

(Explosion-proof class) Ex d IIC T6 X

(Ambient Temperature)* -20°C to +53°C

(Electrical Data) Supply voltage: DC24V 45mA


Detecting element output signal: DC0.1V 0.2mA

Analog signal output: DC24V 22mA

(Applicable Standard) JNIOSH-TR-NO. 43(2008)

*The ambient temperature on explosion-proof certification is described above and it indicates the temperature range which can maintain the explosion-proof performance and not the product performance, see "10. Product specifications" for details.

Installation




3


Product Components

3-1. Main unit and standard accessories

<Main Unit> (including a cable gland or an adapter)
 [Export Specification]



[Japan Specification]



<Standard Accessories>

- Operating manual
- Dedicated handling lever (used for the wiring)
- Dedicated control key




CAUTION

Use the supplied dedicated control key to operate the detector. If products other than these accessories are used, key operations cannot be accepted properly.

Note that the control key is made of extremely strong magnet. Putting it close to a magnetic product, such as a credit card and ID card, may damage the stored data.


3-2. Names and functions for each part



(1)	MENU/ESC key	Used to enter the maintenance mode. It is also used to cancel in a specific mode.
(2)	SET key	It is used for value confirmation and so on in a specific mode.
(3)	▲ key	Used to switch menus or change a value (UP).
(4)	▼ key	Used to switch menus or change a value (DOWN).
(5)	Power lamp	Power lamp. Detection mode: It lights up in green. Maintenance mode: It blinks in green.
(6)	Alarm lamp	Alarm lamp. It lights up in red when the alarm setpoint value is reached.
(7)	Fault lamp	Fault lamp. It lights up in yellow when an abnormality is detected in the detector.
(8)	Concentration value display	Display the gas concentration and so on.

3-3. Block diagram

<Electric Diagram>



4

How to Use

4-1. Before using the gas detector

Not only the first-time users but also the users who have already used the detector must follow the operating precautions.

Ignoring the precautions may damage the gas detector, resulting in inaccurate gas detection.

4-2. Precautions for installation sites



CAUTION

This is a precision device. Because the detector may not provide the specified performance in some places (environments), check the environment in the installation point, and then take appropriate actions if necessary.

Because the detector plays an important role for safety and disaster prevention, as many units of the detector as needed must be installed in appropriate points.

Because points where gases leak and remain easily are different depending on the types of gases and the working areas, please decide carefully on installation points and the number of units to be installed.

Do not install the detector in a place with vibrations or shocks.

The detector consists of sensitive electronic parts. The detector must be installed in a stable place without vibrations or shocks and it cannot drop.

Do not install the detector in a place exposed to water, oil or chemicals.

When selecting installation points, avoid a place where the detector is exposed to water, oil or chemicals.

Do not install the detector in a place where the range of operating temperatures is exceeded.

The detector must be installed in a stable place where the operating temperature is maintained and does not change suddenly.

<Export Specification> -10 to +40°C

<Japan Specification> -10 to +40°C

Do not install the detector in a place exposed to direct sunlight or sudden changes in the temperature. When selecting installation points, avoid a place where it is exposed to direct sunlight or radiant heat (infrared rays emitted from a high-temperature object), and where the temperature changes suddenly. Condensation may be formed inside the detector, or the detector cannot adjust to sudden changes in the temperature.

Keep the detector (and its cables) away from noise source devices.

When selecting installation points, avoid a place where high-frequency/high-voltage devices exist.

Do not install the detector in a place where maintenance of the detector cannot be performed or where handling the detector involves dangers.

Regular maintenance of the detector must be performed.

Do not install the detector in a place where the machinery must be stopped when maintenance is performed in its inside, where parts of the machinery must be removed to perform maintenance, or where the detector cannot be removed because tubes or racks prevent access to it. Do not install the detector in a place where maintenance involves dangers, for example, near a high-voltage cable.

Do not install the detector in machinery which is not properly grounded.

Before installing the detector in machinery, the machinery must be grounded properly.

Do not install the detector in a place where interference gases exist around it.

The detector must not be installed in a place where interference gases exist around it.



WARNING

An oxygen detector is affected by pressure. If the detector is used at a pressure other than atmospheric pressure, a calibration is required under that pressure. The maximum allowable pressure is 110 KPa. The detector cannot be used under reduced pressure. Using the detector at a pressure exceeding the allowable range may cause a trouble in the oxygen sensor performance.

4-3. Precautions for system designing



CAUTION

An unstable power supply and noise may cause malfunctions or false alarms.

The descriptions in this section must be reflected on the designing of a system using the detector.

Using a stable power supply

The external output of the detector may be activated when the power is turned on, when momentary blackout occurs, or while the system is being stabilized. In such cases, use a UPS (uninterruptible power supply), or take appropriate actions on the receiving side.

The detector must be provided with the following power supply.

Power supply voltage	24 VDC±10%: Terminal voltage of the detector			
	Up to 10 milliseconds	Example of actions		
Allowed time of (To recover from the momentary		To ensure continuous operation and		
momentary blackout blackout for 10 milliseconds or more,		activation, install a UPS outside the		
	restart the detector.)	detector.		
	Do not use it with a power supply of	Example of actions		
Others	large power load or high-frequency	Use a line filter to avoid the noise source		
noise.		if necessary.		

Introducing protective measures against lightning

If cables are installed outside the factory/plant, or if internal cables are installed in the same duct as the cables coming from outside the factory/plant, "lightning" will cause problems. Because lightning acts as a large emission source while cables act as a receiving antenna, devices connected to the cables may be damaged.

Lightning cannot be prevented. Cables installed in a metal conduit or under the ground cannot be completely protected from inductive lightning surge caused by lightning. Although complete elimination of disasters caused by lightning is impossible, the following protective measures can be taken.


4 How to Use 4-4. How to install

	Take appropriate measures in accordance with the importance of the facilities and the			
	environment.			
Protection against	Provide protection by a lightning arrester (cable arrester).			
lightning	(Although inductive lightning surge can be transmitted through the cable, it is			
iigiiiiiig	prevented by installing a lightning arrester before the field devices and central			
	processing equipment. For information on how to use a lightning arrester,			
	please contact the manufacturer.)			
Grounding	In addition to lightning, there are more sources of surge noise. To protect units from			
Grounding	these noise sources, the units must be grounded.			

* The lightning arrester has a circuit to remove a surge voltage which damages field devices, so that signals may be attenuated. Before installing a lightning arrester, verify that it works properly.

4-4. How to install

<Installation Dimensions and Maintenance Space>



The following installation requirements must be met to install the detector.

• Attach the detector on the wall and others using four M5 screws.



CAUTION

<u>Do not install the detector in a place where maintenance of the detector cannot be performed or where handling the detector involves dangers.</u>

Regular maintenance of the detector must be performed.

Do not install the detector in a place where the machinery must be stopped when maintenance is performed in its inside, where parts of the machinery must be removed to perform maintenance, or where the detector cannot be removed because tubes or racks prevent access to it. Do not install the detector in a place where maintenance involves dangers, for example, near a high-voltage cable.

4 How to Use 4-4. How to install




WARNING

Install the detector head in the specified direction.
 Installing it in a wrong direction may cause abnormal readings.

• Install the detector while giving due consideration to the conditions for the detection area. Normally, when the detector is used to detect oxygen deficiency, the detector head should be installed at a face level.

- When there is a risk of oxygen deficiency due to a gas heavier than the air (CO2 etc.), install
 the detector in consideration for the specific gravity of the gas.
- Do not install the detector in a damp place where condensation occurs. Water drops formed on the sensor may cause abnormal detection.

Detector head installation direction




4-5. Installation procedure

Attaching External Cable>

[Export Specification]

The export specification comes with adapter A (NPT 1/2) as standard. Please prepare the Ex cable gland that is suitable for the adapter.




[Japan Specification]

Cable gland, Rubber seal, Washer and Eccentric washer are mounted on SD-1OX main body as shown below.

Attach the parts:

Pass a Cable (e.g. CVVS) through the Cable gland, Eccentric washier, Washer, and Rubber seal, and drawn into detector inside. (For wire methods, see "4-6 How to wire".)

Tightens Cable gland, and then tightens hexagon socket set screw (M4×6), and fix cable gland.





WARNING

• If a scratch, crack or other abnormality is found on the container or explosion-proof structure joint surface, stop using the detector immediately and contact RIKEN KEIKI.



CAUTION

- Tighten 40 N·m or larger torque value to ensure the tightening cable gland and seal plug.
- If it is difficult to tighten the cable gland, grease its screw part and then tighten it with tools.

4-6. How to wire



CAUTION

• Be careful not to damage the internal electronic circuit when wiring. In addition, be careful not to apply stresses on the detector when (overweight) cables are installed.


- The power cables and signal cables must not be installed together with the motor power cables
 etc. When these cables must be installed together for unavoidable reasons, put the power
 cables and signal cables in a metal conduit. The conduit must be connected to a grounding
 circuit.
- When stranded wires are used, prevent wires from contacting each other.
- Use the dedicated handling lever to wire.
- Use appropriate cables to wire.

<Recommended Cables>

3-wire	CVVS 1.25mm ² or 2.0mm ² /3-core
4-wire	CVVS 1.25mm ² or 2.0mm ² /4-core

<Parts. length table of outside conductor lead-in>

Cable overall outer diameter(mm)	Rubber seal inner diameter(mm)	Washer inner diameter(mm)
From Φ9.6 to less than Φ11.0	Ф11	Ф12
From Φ11.0 to less than Φ11.5	Ф12	Ф12
From Φ11.5 to less than Φ12.0	Ф12	Ф14
From Φ12.0 to less than Φ13.0	Ф13	Ф14
Ф13.0	Ф13.5	Ф14




NOTE=


• The following table shows an example of overall outer diameters of cables. Use them for reference. The overall outer diameters must be checked because they somewhat vary between manufacturers.

Number of core	CVV 1.25mm ²	CVV 2mm ²	CVVS 1.25mm ²	CVVS 2mm ²
2	ø9.5	ø10.5	ø10.0	ø11.0
3	ø10.0	ø11.0	ø10.5	ø11.5
4	ø10.5	ø11.5	ø11.0	ø12.0
5	ø11.5	ø12.5	ø12.0	ø13.0
6	ø12.5	ø13.5	ø13.0	ø14.0

<Figure of Terminal Plate>



24 VDC	DC+	1
24 VDC	DC-	2
4 - 20 mA	Sig+	3
HART	Sig-	4



NOTE •

The 2p of terminal plate is being deemed to be power and the 4-20mA output common (-) terminal, and using as 3-wire is possible.

<Specifications of Terminal Plate>

Specifications of terminal plate

• Rated voltage: 250 VAC

• Rated current: 12 A

However, it depends on cables to be used.

Connection conditions

• Cables: 0.25 - 2.5 mm²

• Bare wire length: 8 - 9 mm

• Connecting tool: Dedicated handling lever (accessory) or driver (edge 3.5 x 0.5 mm)




CAUTION

The specified bare wire length must be observed when the wire insulation is peeled off.

Improper clamping of the wire due to a shorter bare wire length may cause defective electrification or heating.

Catching the wire insulation due to a shorter bare wire length may cause defective electrification or heating.

Exposing the wire due to a longer bare wire length may cause defective insulation or a short circuit. Be careful not to break up the wire. If the wire is broken up when inserted to the terminal, this may cause defective insulation or heating.



Compatible bar terminal

For a bar terminal, the following items are available.

- Bar terminal (ferrule): Model 216 Series (manufactured by WAGO)
- Crimping tool: Model VarioCrimp 4 (206-204) (manufactured by WAGO)



CAUTION

A bar terminal of the specified model must be used. Using other bar terminals invalidates the warranty of the performance.

<How to Connect to Terminal Plate>

When cables are connected to the connectors, use the dedicated lever or a flathead screwdriver to do it as shown below.



CAUTION

The right tools must be used.

In principle, one wire can be connected to one wiring hole.

When the wire is inserted into the driver slot by mistake, it does not contact the conductive part. This may cause defective electric conduction or heating.


When the wire is inserted under the spring by mistake, it does not contact the conductive part. This may cause defective electric conduction or heating.

NOTE =

<How to Use the Dedicated Handling Lever>



Push the lever with your finger to lower the spring in its inside.



While holding down the lever, insert the wire into the (round) wiring hole until it reaches the deepest point. Once the lever is released, the wire is secured.

To check whether the wire is connected securely, pull the wire gently. (Do not pull the wire strongly.)

<Grounding>

Connect the detector to your grounding terminal with the external terminal







WARNING

- Before turning on the detector, do not forget to connect it to a grounding terminal.
- For stable operation of the detector and safety, it must be connected to a grounding terminal. Do not connect the grounding wire to a gas pipe.
- The grounding must be made as D type grounding (below 100 Ω of grounding resistance).
- For the grounding wire, use cable lugs to safely connect it to a grounding terminal without looseness or twist.
- Use ring terminals to connect the grounding terminal with the ground and use the grounding wire with cross-sectional area of 4 mm2 or more for the external grounding terminal.

<Wiring Example>







If HART communication is used, check so that wiring load resistance that constitute 4-20mA loops, amount to 250 Ω -300 Ω .

5

How to Operate


5-1. Preparation for start-up

Before connecting a power supply, read and understand the following precautions. Ignoring these precautions may cause an electric shock or damage the detector.

- Connect the detector to a grounding circuit.
- Check that the wiring is connected to external device properly.
- Check that the power supply voltage is compliant with the specifications.

5-2. Basic operating procedures

Normally, the detection mode is used for normal operations. (The detection mode is activated after the power is turned on.)





WARNING

When the detector enters other mode from the detection mode while an alarm is activated, the alarm is reset.

5-3. How to start the gas detector

- Before supplying power (24 VDC) to the detector, check that the detector is installed properly.
- Supply power (24 VDC) to the detector (turn on the detector).
- After the detector completes the start-up, it enters the detection mode swiftly.




CAUTION

- Do not turn off the detector during the initial clear. The detector is reading the internal memory during the initial clear.
- If the detector is installed newly or the new sensor is replaced, the sensor must be warmed up
 for a specified period which is determined depending on the type of the sensor after the
 detector is started.
- After the warm-up is completed, perform a span adjustment.

NOTE =

<<Start-up Operation Procedures>> (approximately 25 seconds for system check of the detector and alarm deactivation)

Power on -> Initial clear (approximately 25 seconds) -> Detection mode



5 How to Operate 5-4. Modes

5-4. Modes

Details on each mode are provided as follows.



CAUTION

• Do not change the settings if not necessary. Changing the settings without understanding the specifications may cause malfunctions.

Mode	Item	LED display	Details	
Detection Mode	-	Gas concentr ation	Normal state	
	ROM/SUM Display	1-0	Display the program version and others. This is not typically used by the user.	
User Mode	Span Adjustment	1-1	Perform the span adjustment.	
O3CI WIOGC	Setting Display	1-2	Display various setting values.	
	Switch to Maintenance Mode	1-3	Switch to the regular maintenance mode.	
	Test Mode	2-0	Perform various tests. 2-0.0 Gas Test 2-0.1 Alarm Test 2-0.2 Fault Test 2-0.3 LED Test 2-0.4 Memory Test	
	Zero Adjustment (N2 adjustment)	2-1	Perform the zero adjustment(N2 adjustment)t.	
	Span Adjustment (air adjustment)	2-2	Perform the span adjustment(air adjustment).	
Maintenance Mode (Regular maintenance)	Zero/Span Initialization	2-3	Initialize zero/span values.	
	Environmental Setting	2-4	Used for various environmental settings. 2-4.0 Sensor Power Supply ON/OFF 2-4.1 INHIBIT Setting 2-4.2 Alarm Setpoint Setting 2-4.3 Alarm Delay Time Setting 2-4.4 Alarm Pattern Setting 2-4.5 Suppression Type Setting 2-4.6 Suppression Value Setting 2-4.7 Not used. 2-4.8 Not used. 2-4.9 Not used. 2-4.9 Not used. 2-4.A Maintenance Mode External Output Setting 2-4.B External Output Adjustment 2-4.C Alarm Test External Output Setting 2-4.D Not used. 2-4.E Password Setting 2-4.F Sensor Fault Alarm Pattern Setting 2-4.H Not used.	
	Display	2-5	Display various electrical settings. This is not typically used by the user.	
	Switch to Factory Mode	2-6	Not used.	
	Switch to User Mode	2-7	Returns to the user mode.	

5 How to Operate 5-5. User mode

5-5. User mode



WARNING

After the adjustment is completed, never fail to press MENU/ESC key to return to the detection

(If the detector remains in the user mode, it automatically returns to the detection mode in ten hours.)


Detection Mode

Press the MENU/ESC key for three seconds.





1-0. ROM/SUM Display Display the program version and others. This is not typically used by the user.




1-1. Span Adjustment (air adjustment) Perform the span adjustment.

1-2. Various Setting Display Display various setting values.

1-3. Maintenance Mode Switching

Switch to Regular maintenance

mode.



Setting Display => P26 0 **SET**

SET

3 1-


See "Regular maintenance mode". **SET**

Span Adjustment(air adjustment)

=> P25

5 How to Operate 5-5. User mode

This is used to perform the span adjustment(air adjustment).




NOTE -

• If the span adjustment(air adjustment) failed since the span point was significantly fluctuated from around 20.9%, it returns to 1-1 after FAIL rather than PASS is displayed. In this case, the span adjustment(air adjustment) has not been completed.

5 How to Operate 5-6. How to exit

<Setting Display "1-2">

Display various setting values.



5-6. How to exit

To turn off the detector, turn off the power supply (24 VDC) to the detector.



WARNING

- When the detector is turned off, an alarm may be triggered on the upper (central) system.
- Before turning off the detector, INHIBIT (point skip) on the upper (central) system must be
 activated. Decide whether the power can be turned off by checking the operation of the devices
 connected to the external output or external contact output terminal of the detector.

6

Operations and Functions

6-1. Alarm activation

Oxygen concentration alarm: Activated when the detected oxygen concentration reaches or exceeds the alarm setpoint value. <<Auto-Reset Operation>>

NOTE -

The alarm setpoint is factory-set. Although the alarm delay time (standard: 2 seconds) works in the detector to prevent a false activation, it can be cancelled if not needed.

<Display Operation>

Oxygen Concentration Display

In case of over the detection range (Over Scale), " $\cap \cap \cap \cap$ " is displayed on the LED.

Power Indicator Lamp (POWER: Green)


During operation, this lights up continuously.

Alarm Indicator Lamp (ALM: Red)

This lights up when the alarm setpoint value is reached or exceeded.

<Alarm Activation>

It lights up when the oxygen concentration reaches or exceeds the alarm setpoint value. The alarm activation is reset automatically when the gas concentration drops above the alarm setpoint value.



<Response to Alarm>

An oxygen concentration value exceeds the alarm setpoint

When an alarm is triggered, take actions in accordance with your management rules of alarm. Normally, take the following actions.

Check the reading of the detector.

NOTE •

If a gas leak is momentary, the reading may already have dropped when checking it. In addition, when the alarm is triggered by noise or other incidental conditions other than a gas, the reading may have already dropped.

- Based on your management rules of gas alarm, no one can be allowed to access the monitored zone to ensure safety.
- If the gas concentration display continues to be displayed, close the main valve of the gas, and then
 check that the gas concentration reading dropped.
- Access the gas leak point, equipped with a protective gear to avoid dangers caused by possibly remaining gases, and check whether gases remain or not by using a portable gas detector.
- Check that the point is free from dangers, and take actions to fix the gas leak.

6-2. Fault alarm activation

A fault alarm is triggered when the detector detects abnormalities. After a fault alarm is triggered, the fault lamp (yellow) lights up and an error message is displayed on the LED, determine the causes and take appropriate actions.

After the detector is successfully returned from the fault, it restarts with the process normally performed right after it is turned on (initial clear).

If the detector has problems and is repeatedly malfunctioning, contact RIKEN KEIKI immediately.

<Display Operation>

Fault Detail Display

Display a message indicating the fault detail.

Fault Lamp (FAULT: Yellow)

This lights up when a failure occurs.

NOTE


For information on malfunctions (error messages), see "9. Troubleshooting".

6-3. External output operation

Signa	al Transmission System	Electric current transmission (non-isolated) 4 – 20 mA	
Transmission Path		CVVS	
Trans	smission Distance	CVVS 1.25mm ² : Within 1.25km	
Conn	ection Load Resistance	Below 300 Ω	
1	Detection Mode(No alarm)	4 - 20 mA (concentration output)	
2	Detection Mode(Gas alarm)	4 - 20 mA (concentration output)	
3	Initial Clear	Depending on the setting of the section 4. 2.5 mA setting: 2.5 mA 4 mA, HOLD, 4 - 20 mA setting: Air value(17.4 mA)	
4	Maintenance Mode	2.5 mA setting: 2.5 mA 4 mA setting: Air point(17.4 mA) HOLD setting: The previous value retained 4-20 mA setting: 4 - 20 mA (concentration output)	
5	Alarm Test Mode	Output ON setting: 4 - 20 mA (concentration output) Output OFF setting: The previous value retained	
6	Fault Alarm	0.5 mA (Fixed)	
7	INHIBIT	Depending on the setting of the section 4. 2.5 mA setting: 2.5 mA 4 mA, HOLD, 4 - 20 mA setting: Air value(17.4 mA)	
8 Power Off 0 mA		0 mA	

Example of Gas Concentration and External Output

0 - 25 vol% (Oxygen deficiency alarm specification)



CAUTION

- The 4 20 mA output is already adjusted. In case of over scale, an output will not exceed 22 mA.
- Output during INHIBIT or initial clear is based on 4 20 mA output setting in the maintenance mode.


In particular, when the detector is started or the specification is changed, be careful about 4-20 mA output setting. Understand how the detector functions, and take actions, if necessary, so that the signal receiver side can prepare to avoid false alarms.

6-4. Other functions

<Suppression Function>

The sensors used with the detector are influenced by environmental changes (temperature, humidity and other characteristics) in no small measure, which affects the reading. Therefore, the reading might be fluctuated even in a normal environment.

This function obscures influences by environmental changes or other factors that have no meaning for your management rules of alarm. This function is used to hide (suppress) the fluctuation of the reading under the setting value, indicating 20.9 vol%.



NOTE -

- The suppression function is factory-set. The standard setting value is ±0.5 vol% (for 20.9 vol%).
- In the maintenance mode, this function is disabled and the fluctuation of the reading under the setting value is displayed.



CAUTION

A reading under zero is suppressed with the 10% FS suppression.

A reading that gets 10% FS or more under zero is displayed as "-0.0", which prevents an accurate gas detection and needs the zero adjustment.

7

Maintenance

The gas detector is an important instrument for the purpose of safety.

To maintain the performance of the gas detector and improve the reliability of safety, perform a regular maintenance.

7-1. Maintenance intervals and items

- Daily maintenance: Perform maintenance before beginning to work.
- Monthly maintenance: Perform maintenance on the alarm circuit (alarm test) once a month.
- Regular maintenance: Perform maintenance once or more for every six months to maintain the performance as a safety unit.

Maintenance item	Maintenance content	Daily mainte- nance	Monthly mainte- nance	Regular mainte- nance
Power Supply Check	Check that the power lamp lights up.	0	0	0
Concentration Display Check	Check that the concentration display value is zero(oxygen monitor:20.9vol%). When the reading is incorrect, perform the zero adjustment (or span adjustment) after ensuring that no interference gases exist around it.	0	0	0
Alarm Test	Inspect the alarm circuit by using the alarm test function.	I	0	0
Span Adjustment	Perform the span adjustment by using the calibration gas.	1	-	0
Gas Alarm Check	Check the gas alarm by using the calibration gas.	-	1	0

<About Maintenance Services>

- We provide services on regular maintenance including span adjustment, other adjustments and maintenance.
 - To make the calibration gas, dedicated tools, such as a gas cylinder of the specified concentration and gas sampling bag must be used.
 - Our qualified service engineers have expertise and knowledge on the dedicated tools used for services, along with other products. To maintain the safety operation of the detector, please use our maintenance service.
- The followings are typical maintenance services. For more information, please contact RIKEN KEIKI.

Main Services

Power Supply

Check

Checks the power supply voltage. Verifies that the power lamp lights up.

(Verifies that relevant points can be identified on the system.)

(When a UPS (uninterruptible power system) is used, checks the operation with the UPS.) Verifies that the concentration display value is zero (or 20.9 vol% on the oxygen deficiency

Concentration Display Check

meter) by using the zero gas.

Performs the zero adjustment (fresh air adjustment) if the reading is incorrect.

Flow Rate Check

Checks the flow rate indicator to find abnormalities.

Checks the flow rate by using an external flow meter to verify the correctness of the flow rate indicator on the gas detector. If the flow rate is incorrect, performs the flow rate adjustment.

Filter Check

Checks the dust filter for dust or clogging. Replaces a dirty or clogged dust filter.

Alarm Test

Inspects the alarm circuit by using the alarm test function. Checks the alarm lamps. (Checks ALM1 activation.)

Checks the external alarm. (Checks the activation of the external alarm, such as a

buzzer.)

Span

Check

Performs the span adjustment by using the calibration gas.

Adjustment Gas Alarm

Checks the gas alarm by using the calibration gas.

Checks the alarm. (Checks triggering of alarm when the alarm setpoint is reached.)

Checks the delay time. (Checks time to delay until the alarm is triggered.)

Checks the alarm lamps. (Checks ALM1 activation.)

Checks the external alarm. (Checks the activation of external alarms, such as a buzzer and reset signal.)

Cleaning and Repair of Device

Checks dust or damage on surface, cover or internal parts of the detector, cleans and

repairs such parts of the detector.

Replaces parts which are cracked or damaged.

(Visual Diagnosis)

Device

Uses the keys to check the operation of functions and parameters.

Operation Check

Replacement

of Consumable

Parts

Replaces consumable parts, such as a sensor, filter and pump.

7-2. Maintenance Mode(Regular maintenance)




WARNING


After the adjustment is completed, never fail to press MENU/ESC key to return to the detection mode

(If the detector remains in the maintenance mode, it automatically returns to the detection mode in ten hours.)

Mode	Item	LED display	Details
Maintenance Mode (Regular maintenance)	Test Mode	2-0	Perform various tests. 2-0.0 Gas Test 2-0.1 Alarm Test 2-0.2 Fault Test 2-0.3 LED Test 2-0.4 Memory Test
	Zero Adjustment (N2 adjustment)	2-1	Perform the zero adjustment(N2 adjustment).
	Span Adjustment (air adjustment)	2-2	Perform the span adjustment(air adjustment).
	Zero/Span Initialization	2-3	Initialize zero/span values.
	Environmental Setting	2-4	Used for various environmental settings. 2-4.0 Sensor Power Supply ON/OFF 2-4.1 INHIBIT Setting 2-4.2 Alarm Setpoint Setting 2-4.3 Alarm Delay Time Setting 2-4.4 Alarm Pattern Setting 2-4.5 Suppression Type Setting 2-4.6 Suppression Value Setting 2-4.7 Not used. 2-4.8 Not used. 2-4.9 Not used. 2-4.9 Not used. 2-4.A Maintenance Mode External Output Setting 2-4.B External Output Adjustment 2-4.C Alarm Test External Output Setting 2-4.D Not used. 2-4.E Password Setting 2-4.F Sensor Fault Alarm Pattern Setting 2-4.H Not used.
	Display	2-5	Display various electrical settings. This is not typically used by the user.
	Switch to Factory Mode	2-6	Not used.
	Switch to User Mode	2-7	Returns to the user mode.



<Test Mode "2-0">






WARNING

- Before starting the alarm test (transmission test), provide a notification to the related sections so
 that they can prepare for false abnormalities (external output signals). After the test is completed,
 never fail to press MENU/ESC key to return to the detection mode. (If the detector remains in the
 alarm test mode, it automatically returns to the detection mode in ten hours.)
- In the alarm test mode, the gas concentration equal to the reading is output even in the maintenance mode depending on the setting. Take actions, if necessary, so that the signal receiver side can prepare to avoid false alarms.

<Gas Test "2-0.0">


2-0.0 Press SET key.

Introduce the test gas and perform the gas test.



Perform the gas alarm test by increasing or decreasing the reading with ▲ or ▼ key.

Press MENU/ESC key to cancel the test and to go back to the original state.






WARNING

Because the contact and gas concentration output can be activated by a gas test even in the maintenance mode depending on setting, be careful to perform the test.

<Alarm Test "2-0.1">

2-0.1 Press SET key.

Stop introducing the test gas. When the reading drops, press MENU/ESC key to cancel the test and to go back to the original state.






WARNING


Because the contact and gas concentration output can be activated by an alarm test even in the maintenance mode depending on setting, be careful to perform the test.

<Fault Test "2-0.2">



Fault Test ON/OFF

Select either ON/OFF.
Switch ON and press SET key to activate the fault alarm.
Return to OFF and press SET key to deactivate the fault alarm.
(Pressing MENU/ESC key also enables to cancel this menu and to go back to the original state.)






WARNING

Because the contact and gas concentration output can be activated by a fault test even in the maintenance mode depending on setting, be careful to perform the test.

<LED Test "2-0.3">




LED Test ON/OFF

Select either ON/OFF.
Switch ON and press SET key to turn


on the LED (POWER, ALM, FAULT) and to display the gas concentration of 8.8.8.8...

Return to OFF and press SET key to deactivate this menu.

(Pressing MENU/ESC key also enables to cancel this menu and to go back to the original state.)




<Memory Test "2-0.4">



When CAL. is displayed, the memory diagnosis is performed.

When memory is correct as a result of the diagnosis, PASS is displayed. Press SET key to go back to the original state.



LED

2-0.4

F


 \bigcirc


NOTE =

When abnormalities are found in memory as a result of the diagnosis, PASS is not displayed. F-RA is displayed for the RAM fault, and F-RO is displayed for the ROM fault. For both the RAM fault and ROM fault, F-RA and F-RO are displayed alternately.

<Environmental Setting "2-4">

Set various operations and functions in the environmental setting.






7-3. Gas calibration method

Perform a gas calibration in each mode (zero adjustment mode and span adjustment mode) using the calibration gas.

- Zero adjustment gas(more than N2 99.9vol% collected in a gas sampling bag)
- Span adjustment gas (Air)
- Gas sampling bags



This is used to perform the span adjustment (air adjustment).




WARNING

When the span adjustment (air adjustment) is performed in the atmosphere, check the atmosphere for freshness before beginning the adjustment. If interference gases exist, the adjustment cannot be performed properly, thus leading to dangers when the gas leaks.



CAUTION

Before starting the span adjustment (air adjustment), provide the detector with air and wait until
the reading is stabilized.



NOTE

• If the span adjustment failed since the reading was significantly fluctuated from the introduced gas concentration, it returns to 2-2 after FAIL rather than PASS is displayed. In this case, the span adjustment has not been completed.

Check if the correct concentration calibration gas has been provided and check the gas adjustment tool for a gas leak. If there are mistakes, perform the span adjustment again.

If the span adjustment cannot be performed even when there is no mistake or after recalibration, the gas sensor life might have expired. After turning off the power supply, please contact RIKEN KEIKI.


<Zero Adjustment (N2 adjustment) "2-1">

This is used to perform the zero adjustment(N2 adjustment).

When the detector is used to prevent oxygen deficiency, zero adjustment is not required.

NOTE!

 Before starting the zero adjustment(N2 adjustment), provide the detector with the zero adjustment gas(more than N2 99.9vol%) and wait until the indicator is stabilized.



NOTE


• If the zero adjustment(N2 adjustment) failed since the zero point was significantly fluctuated from around 0.0, it returns to 2-1 after FAIL rather than PASS is displayed. In this case, the zero adjustment(N2 adjustment) has not been completed.

7-4. Parts replacement

<Sensor Replacement>

Gas sensor replacement is according to the following procedures.

- (1) Turn off the power supply (24VDC) to the detector.
- (2) After loosened hexagon socket set screw (2 places) by hex wrench (nominal size 2), remove and turning sensor guard. Make sure that a slip ring is in the sensor guard.
- (3) Remove the flame arrester.
- (4) Remove the sensor, and mount a new sensor.



(5) After replacing a sensor, put on a sensor guard and set hexagon socket set screw (2 places) firmly.



※It is efficient that put a lug in front side and let it.

NOTE -

- Rotate to fit the sensor neatly into detector.
- (6) After sensor replacement, turn on the power and warm up, then do the zero point adjustment and span adjustment(air adjustment) by [7-3. Calibration method].

NOTE

The gas calibration using the standard gas is required after the sensor is replaced. Please request it to our overseas sales department or local representatives.




WARNING

- Do NOT exchange any parts while turn on electricity. It may happen fire.
- Turn off the power supply, and then replace parts.

<Replacement of Regular Replacement Parts>

Replacement of Slip ring, Rubber seal E and O ring, are according to the following procedures.

- (1) After loosened hexagon socket set screw by hex wrench (nominal size 2), rotate the sensor guard, and remove the flame arrester and sensor.
- (2) Replace a Slip ring, Rubber seal E and O ring, and then Sensor, Flame arrester and Sensor guard are sequentially mounted. Then, screw up the Sensor guard by Hexagon socket set screw (2 places).



List of recommended regular replacement parts

No.	Item	Maintenance intervals	Replacement intervals (year)	Quantity (pieces/unit)
1	O ring (for guard)	6 months	3 - 8 years	1
2	O ring (for sensor)	6 months	3 - 8 years	1
3	Rubber seal E	6 months	3 - 8 years	1
4	Slip ring	6 months	3 - 8 years	1

^{*} The operation must be checked after replacement by a qualified service engineer. For the stable operation of the unit and safety, ask a qualified service engineer to take care of replacement of the parts that operation must be checked. Request RIKEN KEIKI for operation check.

NOTE =

• The above replacement intervals are recommendation only. The intervals may change depending on the operating conditions. These intervals do not mean the warranty periods either. The result of the regular maintenance may determine when to replace the parts.

8

Storage, Relocation and Disposal

8-1. Procedures to store the gas detector or leave it for a long time

The gas detector must be stored under the following environmental conditions.

- In a dark place under the normal temperature and humidity away from direct sunlight
- In a place where gases, solvents or vapors are not present

8-2. Procedures to relocate the detector or use it again

When the detector is relocated, select a new place in accordance with "4-2. Precautions for installation sites" and "4-4. How to install".

For information on wiring, see "4-5. How to wire". The unpowered time must be minimized when the detector is relocated.



CAUTION

• When using a relocated or stopped/stored detector again, never fail to perform a calibration. For information on readjustment including a calibration, please contact RIKEN KEIKI.

8-3. Disposal of products

When the detector is disposed of, it must be treated properly as an industrial waste in accordance with the local regulations.

9

Troubleshooting

The Troubleshooting does not explain the causes of all the malfunctions which occur on the detector. This simply helps to find the causes of malfunctions which frequently occur. If the detector shows a symptom which is not explained in this manual, or still has malfunctions even though remedial actions are taken, please contact our overseas sales department or local representatives.

<Abnormalities on Unit>

Symptom/Display	FAULT	Causes	Actions
		The wiring is not correct.	Connect the wiring properly.
		The terminal plate is removed.	Connect the terminal plate properly.
The power cannot be turned on.		Abnormalities/momentary blackout of power supply system	Provide the rated voltage. Take measures such as checking or adding the UPS, power supply line filter and insulation transformer.
		Cable abnormalities (open circuit/not connected/short circuit)	Check the wiring of detector and related devices around it.
Abnormal operations	-	Disturbances by sudden surge noise, etc.	Turn off and restart the indicator/alarm unit. If such a symptom is observed frequently, take appropriate measures to eliminate the noise.
Sensor abnormalities E-1	•	The sensor is not connected or improperly connected.	Check if the cable is securely fastened to the terminal plate.
<u>[1]</u>		Faults of the sensor	Replace the sensor with a new one.
System		The rated voltage is not supplied to the detector.	Check the power supply, and supply the rated voltage.
abnormalities E-9	•	Abnormalities of ROM, RAM or EEPROM inside of the detector	Please contact our overseas sales department or local representatives.

<Abnormalities of Readings>

Symptoms	Causes	Actions
	Drifting of sensor output	Perform the zero adjustment (or span adjustment).
The reading rises (drops) and it remains so.	Slow leak	A very small amount of the gas to be detected may be leaking (slow leak). Because ignoring it may cause dangers, take a remedial measure, i.e., taking actions the same as those for the gas alarm.
	Environmental changes	Perform the zero adjustment (or span adjustment).
A gas alarm is triggered	Disturbance by noise	Turn off and restart the detector. If such a symptom is observed frequently, take appropriate measures to eliminate the noise.
despite of no gas leak and no other abnormalities at the detection point.	Sudden change in the environment	When the environment (temperature etc.) changes suddenly, the detector cannot adjust to it and is affected by it. In some cases, the detector triggers an indication alarm. Because the detector cannot be used under sudden and frequent environmental changes, any preventive actions to eliminate them should be taken by the user.
Slow response	Deteriorated sensor sensitivity	Replace the sensor with a new one.
Span adjustment	Improper calibration gas concentration	Use the proper calibration gas.
impossible	Deteriorated sensor sensitivity	Replace the sensor with a new one.

10

Product Specifications

10-1. List of specifications


[SD-10X(TYPE H)]<ATEX/IECEx/UKEX·Oxygen Deficiency Alarm Specification>

Detection principle	Galvanic cell method
	O2
Detectable gas Gas concentration display	
	LED(4digits • 7segments)
Measuring range	0~25vol%
Detection method	Diffusion method
Alarm preset point	18vol%(L) 【Standard】
Power indication	POWER lamp lighting(green)
Indicate accuracy*1	Within ±0.7vol%
(under an identical condition)	
Response time*1	Within 30sec(T90)
(under an identical condition)	
Alarm accuracy*1	Less than ±1vol%(against alarm preset point)
(under an identical condition)	
Alarm-delay time*1	By anoxia alarm(Alarm preset point:18vol%), Less than 5sec.
(under an identical condition)	(When introducing 10∼11vol% gas)(Without piping delay time and communication
	delay time.)
Gas alarm type	Single alarm(L)
Gas alarm indication	ALM lamp lighting(red)
Gas alarm action	Auto-recover
Trouble alarm · Self diagnosis	System failure(E-9)/Sensor failure(E-1)/Out of operating temperature range(E-27)
Trouble alarm indication	FAULT lamp lighting(yellow)/content display
Trouble alarm action	Auto-recover
Transmission scheme	Three-wire analog transmission(in common with power supply <power common="" signal,="" supply,="">) or Two-wire analog transmission</power>
Specification of	DC4~20mA
Transmission	(linear • load resistance less than 300Ω • resolution ability 200division)
Communication scheme	HART 7
Transmission cable	CVVS worth of shield cable(1.25mm ²)·3-core or 4-core
	CVVS worth of shield cable(2.0mm ²)·3-core or 4-core
Transmission distance	Less than 1.25km in case of CVVS 1.25mm ²
	Less than 2.0km in case of CVVS 2.0mm ²
Various functions	Alarm delay/Suppress
Power supply	DC24V±10%
Power consumption	MAX.1.1W
Cabling port	Flame proof packing method <g3 4=""></g3>
	(Compatible cables ϕ 9.6~10.5mm in outer diameter and rubber seals ϕ 11mm in
	inner diameter <mounted>)</mounted>
	(Compatible cables ϕ 11.5~12.5mm in outer diameter and rubber seals ϕ 13mm in
	inner diameter <accessory)< td=""></accessory)<>
	or Adapter A <npt1 2=""> or Adapter B<npt3 4=""></npt3></npt1>
Initial clear	Approx.25sec
Operating temperature	-10~+40°C(non-rapidly-vary)
Operating humidity	Less than 95%RH(non-condensing)
Structure	Wall mounting type
Explosion-proof structure	Flame proof structure
Explosion-proof grade	II 2G Ex db II C T6 Gb(ATEX/UKEX) / Ex db II C T6 Gb(IECEx)
Outer dimension	Approx.148(W)×208(H)×88(D)mm(projection excluding)
Weight	Approx.2.5kg
vvoigiit	/ Approx.2.ong

Color	Munsell 7.5BG5/2

^{*1} In conformity to JIS T8201 2010(Oxygen deficiency indicator).

Outline Drawings




[SD-1OX(TYPE H)]<Japan Ex•Oxygen Deficiency Alarm Specification>

	Ex-Oxygen Denciency Alarm Specification>
Detection principle	Galvanic cell method
Detectable gas	02
Gas concentration display	LED(4digits • 7segments)
Measuring range	0~25vol%
Detection method	Diffusion method
Alarm preset point	18vol%(L) [Standard]
Power indication	POWER lamp lighting(green)
Indicate accuracy*1	Within ±0.7vol%
(under an identical condition)	VIIIII 20.1 VOI/V
Response time*1	Within 30sec(T90)
(under an identical condition)	
Alarm accuracy*1	Less than ±1vol%(against alarm preset point)
(under an identical condition)	()
Alarm-delay time*1	By anoxia alarm(Alarm preset point:18vol%), Less than 5sec.
(under an identical condition)	(When introducing 10∼11vol% gas)(Without piping delay time and communication
	delay time.)
Gas alarm type	Single alarm(L)
Gas alarm indication	ALM lamp lighting(red)
Gas alarm action	Auto-recover
Trouble alarm · Self diagnosis	System failure/Sensor failure
Trouble alarm indication	FAULT lamp lighting(yellow)/content display
Trouble alarm action	Auto-recover
Transmission scheme	Three-wire analog transmission(in common with power supply <power signal,<="" supply,="" td=""></power>
	common>) or Two-wire analog transmission
Specification of	DC4~20mA
Transmission	(linear • load resistance less than 300Ω • resolution ability 200division)
Communication scheme	HART 7
Transmission cable	CVVS worth of shield cable(1.25mm ²)·3-core or 4-core
	CVVS worth of shield cable(2.0mm²)·3-core or 4-core
Transmission distance	Less than 1.25km in case of CVVS 1.25mm ²
Transmission distance	Less than 2.0km in case of CVVS 2.0mm ²
Various functions	Alarm delay/Suppress/HART communication
Power supply	DC24V±10%
Power consumption	MAX.1.1W
Cabling port	Flame proof packing method <g3 4=""></g3>
Cabining port	(Compatible cables ϕ 9.6~10.5mm in outer diameter and rubber seals ϕ 11mm in
	inner diameter <mounted>)</mounted>
	(Compatible cables ϕ 11.5~12.5mm in outer diameter and rubber seals ϕ 13mm in
	inner diameter <accessory)< td=""></accessory)<>
Initial clear	Approx.25sec
Operating temperature	-10~+40°C(non-rapidly-vary)
Operating lumidity	Less than 95%RH(non-condensing)
Structure	Wall mounting type
Explosion-proof structure	Flame proof structure
Explosion-proof grade	Ex d II C T6 X
Outer dimension	Approx.148(W)×208(H)×88(D)mm(projection excluding)
Weight	Approx.2.5kg
Color	Munsell 7.5BG5/2

^{*1} In conformity to JIS T8201 2010(Oxygen deficiency indicator).

Outline Drawings



NOTE =

For information on other specifications such as leak alarm specifications, please contact RIKEN KEIKI.

10-2. List of accessories

Operating manual · · · · · · one
 Safety manual · · · · · · · one
 Dedicated handling lever · · · one


Dedicated control key · · · · · The supplied quantity depends on the number of units to be delivered.

1 to 10 units	One
11 to 20 units	Two
21 to 50 units	Three
Over 51 units	Four

10-3. Detection principle

[Detection principle]

A negative electrode of noble metal and a positive electrode of lead are placed in a resin container filled with electrolyte. A part of the container is opened and covered with a barrier. The negative electrode is installed in contact with the barrier. Lead wires are drawn from the positive and negative electrodes to obtain outputs. A thermistor is connected between the lead wires to perform temperature compensation for the sensor outputs.



11

Definition of Terms

Galvanic cell type	cell type This is a principle of the sensor installed in the detector head. See "10-3. Detection principle" for details.	
vol% A unit used to express the percentage of a specific substance (or gas) in a volution.		
Atmosphere An atmosphere with a temperature within the range of -10 to 40°C and a humidi 95% RH or less at one atmospheric pressure (1013 hPa).		
Initial clear Output from the detector head fluctuates for a while after turning on the power. This is a function to prevent triggering alarm during that time.		
Full scale Maximum value of the detection range.		
Calibration	Adjusts the readings to the calibration gas concentration value by using the calibration gas.	
Zero suppression A function to cut off the specific drifting that the sensor has.		
Alarm delay time	A function which temporarily suspends activation to prevent a false alarm caused by noise from its outside.	
INHIBIT	The gas detection function is temporarily suspended during maintenance etc. of the detector. This is also called "point skip", which has the same function.	
HART communication	DC4-20mA signal digitals signalling is superimposed, and is methods that transmit multiple signals.	

Revision History

Issue	Revision details	Issue date
0	First issue	12/5/2016
1	Declaration of Conformity	12/23/2019
2	Declaration of Conformity	4/1/2020
3	Safety information, Declaration of Conformity	8/1/2020
4	Safety information, Declaration of Conformity	8/21/2020
5	Sensor Replacement	9/30/2020
6	Declaration of Conformity	11/12/2021
7	UKEX added, CE Declaration of Conformity changed, UKCA Declaration of Conformity added	



EU-Declaration of Conformity



RIKEN KEIKI Co., Ltd. 2-7-6, Azusawa, Itabashi-ku, Tokyo, 174-8744 Japan declare under our sole responsibility that the following product conforms to all the relevant provisions.

Product Name: Smart Transmitter/Gas Detector Head

Model: SD-1OX, SD-1DOX, SD-1EC, SD-1DEC

Council Directives		Applicable Standards
2014/34/EU	ATEX Directive	EN IEC 60079-0:2018 EN 60079-1:2014
2014/30/EU	EMC Directive	EN 50270:2015
2011/65/EU ^[1]	RoHS Directive	EN IEC 63000:2018

^[1]Including substances added by Commission Delegated Directive (EU) 2015/863

EU-Type examination Certificate No.

DEKRA 13ATEX0035 X

Notified Body for ATEX

DEKRA Certification B.V. (NB 0344) Meander 1051,6825 MJ Arnhem P.O.Box 5185,6802 ED Arnhem

The Netherlands

Auditing Organization for ATEX

DNV Product Assurance AS (NB 2460)

Veritasveien 1 1363 Høvik Norway

The marking of the product shall include the following:

II 2 G Ex db IIC T6 Gb

Place: Tokyo, Japan

Date: Jun. 29, 2022 Takakura Toshiyuki

> General manager **Quality Control Center**

I. Talmolice



UK-Declaration of Conformity

Document No. 320UK23009



We, RIKEN KEIKI Co., Ltd. 2-7-6, Azusawa, Itabashi-ku, Tokyo, 174-8744, Japan declare under our sole responsibility that the following product conforms to all the relevant provisions.

Product Name Smart Transmitter/Gas Detector Head Model SD-1OX, SD-1DOX, SD-1EC, SD-1DEC

Regulations	UK designated Standards
Electromagnetic Compatibility Regulations 2016 (S.I. 2016/1091)	BS EN 50270:2015
The Equipment and Protective Systems Intended for Use in Potentially Explosive Atmospheres Regulations 2016 (S.I. 2016/1107) (UKEX)	BS EN IEC 60079-0:2018 BS EN 60079-1:2014
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 (S.I., 2012/3032)	BS EN IEC 63000:2018

UK-Type examination Certificate No.

DEKRA 21UKEX0368X

Approved Body for UKEX

DEKRA Certification UK Ltd (AB8505) Stokenchurch House, Oxford Road,

Stokenchurch, Buckinghamshire HP14 3SX,

United Kingdom

Auditing Organization for UKEX

DNV Business Assurance UK Ltd (AB8501) 4th Floor Vivo Building, 30 Stamford Street,

London SE1 9LQ, United Kingdom

I Talenton

The marking of the product shall include the following

 $\langle \epsilon_x \rangle$

II 2 G Ex db IIC T6 Gb

Place: Tokyo, Japan

Date: Aug. 31, 2023

Takakura Toshiyuki General manager

Quality Control Center